

MD NEEDLE PEN

Our MicroNeedle Pen Can Also Be Used for
Nano Needling



Medical grade microneedling device that can be used with or without the cord. The high quality unit does not come with a steep purchase price or ridiculous disposable cost.

MD NEEDLE PEN

MD Needle Pen Purchase Includes:

- 1 – MD Needle Pen with 2 Lithium batteries and cord.
(Can be used corded or cordless)
- 10 – Replacement Cartridges
- 25 – Patient Brochures (No Nano Needle at this time)
- 1 – Bottle of Pure Hyaluronic Acid
- Operation & Instruction Manual
- Video Tutorial on use of the MD Needle Pen
- Success Portal Access – Treatment Forms, Setting Sheets, Sample Ads (Eblast, FaceBook & Social Media Posts), Before & After Photos, FAQ's
- Adjustable Needle Depth – 0.25mm – 2.5mm
- Adjustable Speed- 3000 – 10,000 rpm
- 1 Year Replacement Warranty
- Discount – On Hands-on Training Class

COST - \$995


[Competitor Comparison](#)



MD NEEDLE PEN

Microneedling Patient Applications:

- Skin Rejuvenation
- Fine Lines & Moderate Wrinkles
- Skin Laxity
- Acne Scars
- Diminished Texture, Tone, Color
- Stretch Marks
- Traumatic Scars
- Photo Aging & Sun Damage
- Hyperpigmentation
- Hypopigmentation
- Hair Regrowth




MD NEEDLE PEN

Replacement Cartridges



Micro Needling

- 12 Needle needle Pins (33 Gauge)
- Surgical Grade Stainless Steel
- Needle depth 0.25 – 2.5 mm
- Anti-back-flow design
- Automated vertical needling to reduce discomfort
- Vented port to reduce pressure and protect skin
- Fluid filter to prevent contamination
- Individually packaged and EO Gas Sterilized



Nano Needling

- 81 mono-crystalline silicone pyramids
- Needle length 0.15 mm
- 1/3 Diameter of human hair
- Semi-conductor fabrication
- Transdermal delivery system
- Treats within the upper epidermis
- No State Board Concerns
- Individually packaged and EO Gas Sterilized

Package of 10 Individual Wrapped Cartridges - \$90

Nano Needling Patient Applications:

- Boost Product Absorption
- Awaken Collagen Producing Cells
- Improve Texture, Tone, Color
- Skin Rejuvenation
- Reduce Fine Lines & Wrinkles
- Diminish Pore Size
- Acne and Acne Scars
- Photo Aging & Sun Damage
- Hyperpigmentation
- Restore Healthy Glow
- Improve Plumpness

Nano Needling procedure uses mono-crystalline silicone pyramids to create nano channels that will enhance the transdermal delivery of active ingredients to effectively treat aesthetic skin conditions.

Create a custom treatment protocol using your existing skin care line and see dramatic change.

This treatment is safe, gentle to the skin, painless, with zero down time. Nano needling is non-invasive and only enters the epidermis and therefore falls into an estheticians scope of practice.

There has been **NO State Board Issues** with this procedure to date.

MD NEEDLE PEN

Training and Certification



Microneedling Class Includes

- Skin physiology
- Patient Selection
- Indications for microneedling use
- Treatment Protocols/Guide lines
- Contraindications - Pre/Post Care
- Skin Preparations
- Disease and infection control
- Safety procedures (needles, blood etc.)
- How to introduce microneedling to your clients
- Advanced Combination treatments
- Hands-on training (procedure technique)
- Certification

8 Hour Class – \$595 or \$495 w/pen purchase
Call for Schedule

MD NEEDLE PEN

Your Success Is Our Success

Educate and Support the Customer

- Provide a quality medical grade microneedling device.
- Keep customer startup and disposable costs low.
- Educate customer on proper use for the most effective results.
- Support the customer (phone, email, video, marketing...)

The mission of MD Needle Pen is to make every customer successful by giving the tools needed to do it. Let us show you how to add this profitable procedure that has one of the highest patient satisfaction rates to your practice.